



Invacare **Standing Transfer Vest**

The **Standing Transfer Vest** is applied around the trunk. The client's feet must be on the floor and the client must also have some standing ability. The **Standing Transfer Vest** is used for those with reduced balance and/or walking ability. The **Standing Transfer Vest** is especially useful for toileting. It also reduces the risk of injury to the carer. A groin support is available as an accessory to the **Standing Transfer Vest**. This helps keep the vest in place.

Invacare **Standing Transfer Vest**

Size	A	B	Product number / Material	
			Polyester fabric	Groin Band
XS	150	750	1485404	–
S	150	850	1485405	1492662
M	180	950	1485406	–
L	180	1050	1484683	1492663
XL	180	1200	1484684	–